



Date: - 18th December, 2018

Corrigendum
For
Automated slide stainer with parallel batches for the
Department of Pathology

NIT Issue Date	: 30 th October, 2018
NIT No.	: Admn/Tender/123/2018-AIIMS.JDH
Pre-Bid Meeting	: 13 th November, 2018 at 03:30 PM
Earlier Last Date of Submission	: 13 th December, 2018 at 03:00 PM
Extended Last Date of Submission	: 27 th December, 2018 at 03:00 PM
Bid opening	: 28 th December, 2018 at 03:15 P.M

The following revised and additional specification will be added:-

- 1. Page No. 10 Heading: Automated slide stainer with parallel batches “3rd Point” :**
For
Equipment should have solvent resistant color touch screen to monitor the staining process using the graphical process representation.
Read
Equipment should have solvent resistant color touch screen control panel with keyboard to monitor the staining process.
- 2. Page No. 10 Heading: Automated slide stainer with parallel batches “5th Point” :**
For
The equipment should have 34 reagent stations and 6 wash stations of 450ml capacity.
Read
The equipment should have 26-34 reagent stations and 6 wash stations of 450ml capacity.
- 3. Page No. 10 Heading: Automated slide stainer with parallel batches “6th Point” :**
For
The equipment should be programmable for 50 programs of upto 40 steps each with incubation time setting from 0sec to 59 minutes 59 seconds.
Read
The equipment should be programmable for 50 programs in multiple steps with incubation time setting from 0 sec to 59 minutes 59 seconds.

- 4. Page No. 10 Heading: Automated slide stainer with parallel batches “7th Point” :**
For
Optional integrated oven with temperature setting from 40°C to 70°C for optimal slide drying is preferred.
Read
Optional integrated oven with temperature setting from 40°C to 65°C for optimal slide drying is preferred.
- 5. Page No. 10 Heading: Automated slide stainer with parallel batches “10th Point” :**
For
Agitation programmable from 0-20 times or continuous should be available.
Read
Agitation facility 0-20 times should be available.
- 6. Page No. 10 Heading: Automated slide stainer with parallel batches “11th Point” :**
For
Reagent management system, station information on touch screen and data logging should be available.
Read
Reagent management system, station information on touch screen or control panel and data logging should be available.