

Dated: 04-May-23

CORRIGENDUM

Rate Contract for Procurement of Consumables for Central Sterile Supply Department (CSSD)

NIT No.	:	Admin/RC/09/2022-AIIMS.JDH
NIT Issue Date	:	08-Apr-2023
Earlier Last Date of Submission	:	08-May-2023 up to 03:00 PM
Revised Last Date of Submission	:	23-May-2023 up to 03:00 PM
Revised Date on Opening of Bids	:	24-May-2023 up to 03:00 PM

Sr.	Sr. No.	For	to be Read As
01	Technical Bid 12.I.(VIII) (Page No. 03)	Financial Status: - Bidder must have an average annual turnover of ₹ 10 Crores (Ten One Crores) during the last 3 financial years (2019-20, 2020-21 & 2021-22) for similar products (Copies of Profit & Loss Account and Balance Sheets of Last Three Financial Years (2019-20, 2020-21 & 2021-22) duly authenticate by a Chartered Accountant should be enclosed).	AND Bidder must have an Average Annual Turnover



02	Sr. No. Technical Bid 12.I.(XV) (Page No. 04)	The bidder must upload required Quality Assurance Certifications (If applicable for any item) /documents in techno- commercial bid for each items along with item no. i.e., self- attested copies of CE/USFDA/ISO/BSP/USP/WHO/GMP/BIS as mentioned in the specification of particular tender items, failing which the offer for such items will be rejected.	The bidder must upload required necessary applicable Quality Assurance Certifications for all the quoted items by that bidder in techno- commercial bid. Bidder should upload Quality Assurance Certifications for each item along with item no. i.e., self- attested copies of CE / USFDA / ISO / BSP / USP / WHO / GMP / BIS, failing which the offer for such items will be rejected.
03	32. EMD/BID SECURITY (General Terms & Conditions) (Page No. 12)	EMD / BID SECURITY Annexure	Bank Guarantee Annexure for EMD / Bid Security is attached at the bottom of this document.

	ANNEXURE-XIII		
Sr.	Item Sr. No. (as per Annexure-XIII)	For	to be Read As
04	ltem Sr. No. 01	1 Roll contain 750 nos. of pieces.	The Pack size is of 400-800 labels. AND Quotation should be price per Label basis.
05	ltem Sr. No. 03	Chemical indicator according to EN ISO 11140 – 1 Type 4 indicator to check EO gas penetration. Should be lead free. Should be with adhesive for easy documentation.	Chemical indicator according to EN ISO 11140 – 1 Type 4 indicator to check EO gas penetration. Should be lead free.
06	ltem Sr. No. 04	Self contained BI for EO. EO sterilization processes according to EN ISO 11138- 2,Bacillus Atrophaeus available with population of 106. Compatible to pcd consist of a 4.55 meter long stainless steel tube of 3mm diameter.	Self-contained BI for EO. EO sterilization processes according to EN ISO 11138- 2,Bacillus Atrophaeus available with population of 106. Compatible to re-usable PCD consist of 4.5-5 mtr. Length with tube diameter 2mm-3mm.



07	ltem Sr. No. 05	Process Challenge Device (PCD) to monitor Ethylene Oxide Process, type test according to EN 1422.The PCD consists of a 4.55 m long stainless steel tube of 3 mm diameter connected with a capsule holding the indicator strip sealing one end. As per the requirement of EN 1422	Process Challenge Device (PCD) to monitor Ethylene Oxide Process, type test according to EN1422.The PCD consists of a 4.55 m long stainless steel tube of 3mm diameter connected with a capsule holding the indicator strip sealing one end. As per the requirement of EN1422. Compatible to RE-usable PCD consist of 4.5-5 mtr. Length with tube diameter 2mm-3mm.
08	ltem Sr. No. 06	Process Challenge Device (PCD) to monitor Ethylene Oxide Process, type test according to EN 1422.The PCD consists of a 4.55 m long stainless steel tube of 3 mm diameter connected with a capsule holding the self-contained biological indicator (SCBI) sealing one end. As per the requirement of EN 1422.	Process Challenge Device (PCD) to monitor Ethylene Oxide Process, type test according to EN 1422.
09	ltem Sr. No. 07	1 Roll contain 750 nos. of pieces.	1 Roll should contain 400-800 pieces. AND Quotation should be price per Label basis.
10	ltem Sr. No. 10	Type 4 chemical indicator for steam with adhesive and should full fill requirement of EN ISO-11140-1 and should not change colour in 140°C, 30 Min in dry heat.	Type 4 chemical indicator for steam. Should fulfil requirement of EN ISO-11140-1 and can be used in all sterilization cycle.
11	Item Sr. No. 11	Type 5 chemical indicator for steam with adhesive and should full fill requirement of EN ISO-11140-1 and should not change colour in 140°C, 30 Min in dry heat.	Type 5 chemical indicator for steam Should fulfil requirement of EN ISO-11140-1 and can be used in all sterilization cycles.
12	Item Sr. No. 13	Self contained Biological ampoule should contain bacillus stereothermophilus with the population of 10^5 CFU and type 5 chemical indicator inside the ampoule for instant release should comply to standard EN ISO- 11138-1+3 and type 5 to EN ISO- 11140-1	Self contained Biological ampoule should contain bacillus stereothermophilus with the population of 10^5 CFU and type 5 chemical indicator inside the ampoule for instant release.
13	ltem Sr. No. 14	1 Roll contain 750 nos. of pieces	1 Roll should contain 400-800 pieces. AND Quotation should be price per Label basis.



14	Item Sr. No. 16	indicator according to EN ISO 11140 – 1 Type 4 Indicator to check hydrogen peroxide penetration. Should be lead free with adhesive for easy documentation should monitor 2 min define time and concentration	Indicator according to EN ISO 11140 – 1 Type 4 Indicator to check hydrogen peroxide penetration. Should be lead free. Should monitor 2 min define time and concentration. It should Show colour change. It should be Lead-free.
15	Item Sr. No. 17	Self-contained Biological indicator for hydrogen peroxide sterilization should have Geob.Stearothermophilus with population 10^6 CFU on glass fibre carrier.	Self-contained Biological Indicator /PCD for Hydrogen Peroxide Sterilization. Should have Geo. B. Stearothermophilus with population 10^6 CFU on glass fibre carrier."
16	ltem Sr. No. 21	Compact Process Challenge Device (PCD) for Batch Monitoring System (BMS) - Hollow Load Test according to EN 867 - 5. Type test for large sterilizer as per EN 285 and class B sterilizer as per EN 13060.	Process challenge device (PCD) that contains a 3- hour Readout Biological Indicator. Class 5 Integrator Moving front Chemical integrators to verify all perimeters of steam sterilization together i.e. Time, Temperature and Quality of steam. Designed and certified as per ISO 11138 and EN 866.
17	Item Sr. No. 22	Compact Process Challenge Device - EU version simulating 7 kg cotton pack according to EN ISO 11140-4 and hollow load test according to EN 867-5.	Process challenge device (PCD) that contains a 3- hour Readout Biological Indicator. Class 5 Integrator Moving front Chemical integrators to verify all perimeters of steam sterilization together i.e. Time, Temperature and Quality of steam. Designed and certified as per ISO 11138 and EN 866.
18	ltem Sr. No. 23	Low foam multi enzyme solution. Mild alkaline. PH = 9-10.5 Should work in hard, soft and cold water. Biodegradable, Eco friendly No Eco toxicity CE certified.	Low foam multi enzyme solution. Mild alkaline. PH = 9-10.5 Should work in hard, soft and cold water. Biodegradable, Eco friendly No Eco toxicity CE certified. AND Price to be quoted per Litre basis.



19	ltem Sr. No. 24	4 Litres	The Rate is to be quoted per Litre basis
20	ltem Sr. No. 26	Crepe paper with natural rubber adhesive. For medical use only. Should sustain temperature of 1340c in a autoclave. Should not leave the adhesive behind after sterilization.	Crepe paper with natural rubber adhesive. For medical use only. Should sustain temperature of 134°C or more in an autoclave. Should not leave the adhesive behind after sterilization.
21	ltem Sr. No. 27	Crepe paper with natural rubber adhesive. For medical use only. Should sustain temperature of 1340c in a autoclave. Should not leave the adhesive behind after sterilization.	Crepe paper with natural rubber adhesive. For medical use only. Should sustain temperature of 134°C or more in an autoclave. Should not leave the adhesive behind after sterilization.
22	ltem Sr. No. 30	Complains with the	Complies with the
23	ltem Sr. No. 31	Complains with the	Complies with the
24	ltem Sr. No. 32	Complains with the	Complies with the
25	ltem Sr. No. 33	Complains with the	Complies with the
26	ltem Sr. No. 34	Complains with the	Complies with the
27	ltem Sr. No. 35	Complains with the	Complies with the
28	Item Sr. No. 36	Complains with the	Complies with the



	1		· · · · · · · · · · · · · · · · · · ·
29	ltem Sr. No. 37	Complains with the	Complies with the
30	ltem Sr. No. 38	Addition in specification	ADDED Specification : Product should comply with CDSCO guidelines.
31	Item Sr. No. 39	Addition in specification	ADDED Specification: Product should comply with CDSCO guidelines.
32	ltem Sr. No. 40	Addition in specification	ADDED Specification : Product should comply with CDSCO guidelines.
33	Item Sr. No. 41	Addition in specification	ADDED Specification : Product should comply with CDSCO guidelines.
34	Item Sr. No. 42	Addition in specification	ADDED Specification : Product should comply with CDSCO guidelines.
35	Item Sr. No. 43	Addition in specification	ADDED Specification : Product should comply with CDSCO guidelines.
36	Item Sr. No. 44	Addition in specification	ADDED Specification: Product should comply with CDSCO guidelines.
37	Item Sr. No. 45	60 GSM Paper	60 GSM or More
38	Item Sr. No. 46	60 GSM Paper	60 GSM or More
39	Item Sr. No. 47	60 GSM Paper	60 GSM or More
40	ltem Sr. No. 48	60 GSM Paper	60 GSM or More
41	ltem Sr. No. 49	60 GSM Paper	60 GSM or More
42	Item Sr. No. 50	60 GSM Paper	60 GSM or More
43	ltem Sr. No. 51	60 GSM Paper	60 GSM or More
44	ltem Sr. No. 52	60 GSM Paper	60 GSM or More



45	ltem Sr. No. 69	1 Ltr	Price to be quoted per Litre basis.
46	ltem Sr. No. 70	5 Ltr	Price to be quoted per Litre basis. AND US EPA / US FDA Approved / European Certificates.
47	Item Sr. No. 77	• It should be Spore Strip, self-contained vial • Made up of large number of resistant bacterial spores. • It should be Malfunctioning sterilizers can be repaired and put back into service quickly. • Provide Auto reader free of cost as per requirement of AIIMS Jodhpur. Box of 100 Tests	 * It should be Spore Strip, self-contained vial. Made up of large number of resistant bacterial spores. * It should be Malfunctioning sterilizers can be repaired and put back into service quickly. * Provide Auto reader free of cost as per requirement of AIIMS Jodhpur. * Price to be quoted per Strip / Biological Indicators basis.
48	ltem Sr. No. 79	1) Size 7.5 cm x 7.5 cm. 2) With X-Ray detectable thread. 3) Sterile.	 * Size 7.5 cm x 7.5 cm. * With X-Ray detectable thread. * Sterile * 12 Ply or more
49	ltem Sr. No. 80	1) Size 10 cm x 10 cm. 2) With X-Ray detectable thread. 3) Sterile.	 * Size 10 cm x 10 cm. * With X-Ray detectable thread. * Sterile * 12 Ply or more
50	Item Sr. No. 81	1) Size 7.5 cm x 7.5 cm. 2) Sterile.	 * Size 7.5 cm x 7.5 cm. * Sterile. * 12 Ply or more
51	ltem Sr. No. 82	1) Size 10 cm x 10 cm. 2) Sterile.	 * Size 10 cm x 10 cm. * Sterile. * 12 Ply or more
52	ltem Sr. No. 84	ISO 45001:2015	ISO 14001:2015 And Compatible stainless steel / basket inside the container should be provided.



53	ltem Sr. No. 85	ISO 45001:2015	ISO 14001:2015 And Compatible stainless steel / basket inside the container should be provided.
54	ltem Sr. No. 86	ISO 45001:2015	ISO 14001:2015 And Compatible stainless steel / basket inside the container should be provided.
55	ltem Sr. No. 89	Anionic multienzyme solution with sodium tetraborate decahydrate, should be Biodegradable and free from ethylene glycol. Effective in 3 minutes. Pack size of 3-4L.	Anionic multienzyme solution with sodium tetraborate decahydrate, should be Biodegradable and free from ethylene glycol. Effective in 3 minutes. AND Price to be quoted per Litre basis.
56	ltem Sr. No. 90	Solution containing Benzotriazole, N- (hydroxyethyl) - ethylenediaminetriacetic acid (HEDTA),Dipotassium hydrogen phosphate, Potassium Dihydrogen Phosphate and 0.55% w/v Ortho- Phthalaldehyde. Effective against Human Corona Virus and CytoMegalo Virus. Should also have test strip which consists of sodium sulfite and dyes impregnated and dried on filter paper . PACK SIZE= 5 Litre	Solution containing Benzotriazole, N- (hydroxyethyl) -ethylenediaminetriacetic acid (HEDTA),Dipotassium hydrogen phosphate, Potassium Dihydrogen Phosphate and 0.55% w/v Ortho-Phthalaldehyde. Effective against Human Corona Virus and CytoMegalo Virus. Should also have test strip which consists of sodium sulfite and dyes impregnated and dried on filter paper. AND Price to be quoted per Litre basis.
57	ltem Sr. No. 91	Effervescent 50% troclosene Sodium tablet 5g each tablet.	Effervescent 50% troclosene Sodium tablet 5g each tablet. AND Rate to be quoted per Tablet basis (5G each).



58	ltem Sr. No. 92	100 gm of the granules contains active ingredients more than 40 gm Sodium percarbonate and more than 20 gm tetraacetylethylenediamine and enzymes.	 100 gm of the granules contains active ingredients more than 40 gm Sodium percarbonate and more than 20 gm tetraacetylethylenediamine and enzymes. AND Rate to be quoted per 500 Gram basis.
59	ltem Sr. No. 94	70% v/v Isopropyl alcohol with surfactant and corrosion inhibitors with Hand Sprayer for Rapid surface disinfection. Passes EN 14476, EN13727 and EN13624	70% v/v Isopropyl alcohol with surfactant and corrosion inhibitors with Hand Sprayer for Rapid surface disinfection. Passes EN 14476, EN13727 and EN13624 AND Rate to be quoted per 100 ML basis.
60	ltem Sr. No. 95	Test strips to test MEC level of 2.45 % Glutaraldehyde Solution.	Test strips to test MEC level of 0.55% Orthothalaldehyde Solution. AND Rate to be quoted per Strip basis.
61	ltem Sr. No. 96	Test strips to test MEC level of 0.55% Orthothalaldehyde Solution	Test strips to test MEC level of 2.45 % Glutaraldehyde Solution. AND Rate to be quoted per Strip basis.
62	ltem Sr. No. 109	*Heat sensitive chemical indicator self adherent tape without leaving any adhesive residue, for Steam Sterilization cycles. *Color change should be from White to dark brown for easy identification after processed in Steam Autoclaves of different time and temperature requirements. *Should not change color in dry heat.	 * Heat sensitive chemical indicator self- adherent tape without leaving any adhesive residue, for Steam Sterilization cycles. * Show Colour change for easy identification after processed in Steam Autoclaves of different time and temperature requirements. * Should not change colour in dry heat.



BANK GUARANTEE FORM FOR BID SECURITY

Whereas
against the purchaser's ATE No.
having our registered office at
(Hereinafter called the "Bank") are bound unto AIIMS, Jodhpur (hereinafter called the "Purchaser")
in the sum of for which payment will and truly to be made to the said Purchaser, the Bank binds itself, its successors and assigns by these presents.
Sealed with the Common Seal of thesaid Bank this day of

The conditions of this obligation are:

- 1) If the Bidder withdraws or amends, impairs or derogates from the bid in any respect within the period of validity of this Bid.
- 2) If the Bidder having been notified of the acceptance of his Bid by the Purchaser during the period of its validity:
 - a. If the bidder fails or refuses to furnish the performance security for the due performance of the RateContract / Purchase Orders or
 - b. If the bidder fails or refuses to accept / execute the Contract / Purchase orders or
 - c. If it comes to notice at any time, that the information / documents furnished in its Bid are false orincorrect or misleading or forged

We undertake to pay the Purchaser up to the above amount upon receipt of its first written demand, without the Purchaser having to substantiate its demand, provided that in its demand the Purchaser will note that the amountclaimed by it is due to it owing to the occurrence of one or more the three conditions, specifying the occurred condition(s).

This guarantee will remain in force upto (insert date of additional **Sixty Days After Bid Validity**) and any demand in respect thereof should reach the Bank not later than the above date.

.....

(Signature with date of the authorized officer of the Bank)

.....

(Name and designation of the Officer)

.....

(Seal, name & address of the Bank and address of the Branch)